



PRODUCT DESCRIPTION

A bench process timer with count-up or count-down modes. It is ideal for laboratories, education and industry. Clearly marked non-slip rubber function buttons and a large high contrast LCD display make it easy to use.

- Hours, minutes or seconds
- 12hr or 24hr selectable mode
- 11mm LCD display
- Slide switch to select modes – clock, count-up, timer (for count-down mode)
- Count-down mode with sonic alarm on expiry
- Count-up mode until stopped
- Memory for count-down times
- Impact-resistant ABS case with non slip rubber feet
- Supplied with battery

TECHNICAL SPECIFICATIONS	
Count-Up Capacity	23 hrs, 59 mins, and 59 secs
Count-Down Capacity	23 hrs, 59 mins, and 59 secs with alarm
Resolution	1 second
Accuracy	±1 second in 24 hours
Power	1 x 1.5V AA alkaline battery (LR-6)
Dimensions	85 x 85 x 45mm
Weights	104g