

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

FOUR WAY TIMER & CLOCK


Four timers in one compact unit. 99 hour countdown and up. 12/24 hour clock. Counter automatically continues to count up when countdown time reached.

Product features:

- Four timers in one compact unit.
- Count up and down
- Large, easy to read display.
- Complete with hanging hook, magnetic back and stand.

Product specifications:

Timer keeping:	Hour, minutes, seconds, AM/PM timer indicators (T1 to T4)
Timer count:	99 hours, 99 minutes 12/24 hour clock
Time accuracy standard:	+/- 30 seconds/month
Operating conditions:	0 to 40°C < 95%RH (non-condensing)
Battery:	2 x AAA (included)
CE & RoHS/WEEE:	Compliant
Hazard information (SDS):	See brannan.co.uk for information
LCD dimensions:	23mm x 55mm
Product dimensions:	Net: 70mm x 88mm x 20mm Gross: 110mm x 77mm x 25mm
Product weight:	Net: 86g (including battery) Display box: 17g Gross: 102 .5g



Description	Barcode	Product No
Four way timer and clock	5 011405 282200	28/220/0

