


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
L-Valine - reagent grade, ≥98% (HPLC)

Product Number: **V0500**
CAS Number: 72-18-4
MDL: MFCD00064220
Formula: C₅H₁₁NO₂
Formula Weight: 117.15 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/ml, H ₂ O	
Infrared spectrum	Conforms to Structure
Carbon	50.3 - 52.3 %
Nitrogen	11.7 - 12.3 %
Specific Rotation	26.6 - 28.4 °
(c = 8 in 6 M HCl at 25 deg C)	
Purity (HPLC)	≥ 98 %
Recommended Retest Period	-----
3 years	

Specification: PRD.1.ZQ5.10000033234

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.