



Labworldsoft® 6 Visc

/// Data Sheet

Labworldsoft® 6 Visc software enables the automation of laboratory processes by connecting ROTAVISC viscometers, STARVISC 200-2.5 control or LR-2.ST Starvisc 200 to other devices such as a thermostat or pH Meter through a PC. Measurement processes and parameters can thus be easily recorded for analysis and documentation. In addition, users can e.g. monitor and control the sample temperature and pH measurement.

labworldsoft® 6 Visc at a glance

- > Application for viscosity measurement
- > For all IKA ROTAVISC viscometers, STARVISC 200-2.5 control and LR-2.ST Starvisc 200



designed for scientists

- > For thermostats and pH measuring equipment from numerous manufacturers
 - If your device is missing from the labworldsoft® library, please contact us for possible integration
- > Supports labworldsoft® server
- > Supports the requirements of the FDA 21 CFR Part 11

Hardware and software requirements:

- > Operating System: Windows 10 May 2019 Update (version 1903) 64-bit or later / Windows 11 64-bit
- > .NET Framework: .NET Framework 4.8 or later
- > at least 8 GB RAM
- > at least 10 GB free hard disk space

IKA

designed for scientists