SIGMA-ALDRICH®

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Water - Molecular Biology Reagent

Product Number: CAS Number: MDL: Formula: Formula Weight: W4502 7732-18-5 MFCD00011332 H2O 18.02 g/mol



TEST

Appearance (Turbidity) Appearance (Color)

Appearance (Form) DNAse, Exonuclease Detection NICKase,Endonuclease Detection RNAse Detection Protease Detection by FITC Casein

Specification: PRD.1.ZQ5.10000054299

Specification

Clear Colorless

Liquid None Detected None Detected None Detected

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.