



## PRESS RELEASE

### **New 96-Channel Handheld Electronic Pipette**

The **INTEGRA VIAFLO 96** is a new **96-channel electronic pipette** that enables fast, precise and easy simultaneous transfer of 96 samples.

The VIAFLO 96 has been designed to be used just like a standard handheld pipette. Therefore, no special skills or training are required to operate it. The VIAFLO 96 features the company's popular Touch Wheel™ pipette user interface enabling simple and fast programming of a wide range of pipetting modes including repeat dispense, serial dilute and sample dilute.

Customers will immediately benefit from the increased productivity delivered by their VIAFLO 96. Fast replication or reformatting of 96 and 384 well plates and high precision transferring of reagents, compounds and solutions to or from microplates is now as easy as pipetting with a standard electronic pipette into a single tube.

Four pipetting heads with pipetting volumes up to 12.5 µl, 125 µl, 300 µl or 1250 µl are available for the VIAFLO 96. These pipetting heads are interchangeable within seconds enabling optimal matching of the available volume range to the application performed.

The VIAFLO 96 uses VIAFLO GripTips which provide a perfect seal every time and won't fall off during use. These tips can also be used for other VIAFLO electronic pipettes.

The unique combination of features provided by the VIAFLO 96 will not only increase your productivity but also facilitate accurate delivery of liquids to target vessels, reduce cross contamination, lower operating fatigue and provide consistent results.

For further information on the new VIAFLO 96 please contact Scientific Laboratory Supplies Ltd on telephone 0115 982 1111 or by email [slsinfo@scientific-labs.com](mailto:slsinfo@scientific-labs.com).